


Derma Roller instructions

- Make certain face is clean and dry.
- Use roller gently at first until you know what your skin can tolerate.
- Follow instructions on box for proper direction to roll the AcuLift™ roller.
- Do not use on lips or inside of eye socket (around eyelid).
- Use every 2–3 days until your skin becomes accustomed to the roller.
- If desired, rinse with warm soapy water after use, let dry before closing case.
- Never share with other people, do not use on children or pets.
- Do not use on sunburn, open cuts, active herpes, pustule acne lesions or any other acute infection or inflammation of the skin, have a history of poor wound healing, collagen diseases, blood problems, are pregnant or lactating, are prone to keloid scarring or have diabetes or sore areas.
- Aspirin, Ibuprofen, vitamin E and blood thinning drugs all cause increased bruising and should not be taken in the two weeks leading up to treatment.
- Store in the clear case when not in use.
- Once a week sanitize with Aculift Sanitizing Tablets.



Roll 4 times in each direction.

Possible Side Effects:

Your skin may redden after use, this should go away within a few hours.

You may experience some drying of the skin after use and some flaking. If this happens, stop using the roller until the skin returns to normal (a couple of days).

Use of roller on open wounds may cause infection.

In some individuals, slight bleeding may occur, if this should happen, please area with warm soapy water, apply antiseptic ointment and allow area to heal before continuing use

Use entirely at your own risk. To the maximum extent permitted by law, we are not responsible for negative or damaging results obtained by proper or improper use of this product. Even though many of our customers have experienced very positive results, we cannot guarantee results.